

**Multi-Purpose Benchtop Centrifuges
(Ideal for PRP, Urine or Blood tubes)**



Rotor Types;

Model		FC5000-1B		FC5000-4		T100-A/1		T100-A/2		T100-C	
Rotor Type	Capacity	RPM	RCF(xg)	RPM	RCF(xg)	RPM	RCF(xg)	RPM	RCF(xg)	RPM	RCF(xg)
Angle Rotor	12 x 0.5ml							16000	18320	16000	18320
	48 x 0.5ml									13500	18320
	12 x 1.5ml							16000	18320	16000	18320
	24 x 1.5ml									13500	18320
	4 x 5ml					10000	11000				
	12 x 5ml									15000	18320
	10 x 7ml									13500	18320
	4 x 10/15ml					10000	11000				
	24 x 10/15ml					10000	11000				
	6 x 50ml					10000	11000				
	4 x 100ml					10000	11000				
Swing Rotor	24 x 5ml	4500	4360								
	40 x 5ml	4500	4360								
	24 x 7ml	4500	4360								
	60 x 7ml	4500	4360								
	120 x 7ml			3500	4850						
	88 x 5/7ml			3500	4850						
	8 x 10/15ml	5000	4360								
	12 x 10/15ml	5000	4360								
	16 x 10/15ml	5000	4360								
	24 x 10/15ml	5000	4360								
	28 x 10/15ml			5000	5180						
	32 x 10/15ml	5000	4360								
	64 x 10/15ml			4000	3400						
	72 x 10/15ml			3500	4850						
	32 x 15ml			3500	3400						
	36 x 20ml			3500	4850						
	8 x 50ml	5000	4360								
	16 x 50ml			3500	3400						
	4 x 100ml	5000	4360								
	8 x 100ml			4000	3400						
	4 x 500ml			3500	4850						
Microplate Rotor	2 x 2 x 96 x 0.2ml	5000	4360								
	4 x 2 x 96 x 0.2ml			5000	5180						