

**Yüksek Kapasite ve Yüksek Hızlı Kan
Tüpleri ve PRP için İdeal Santrifüjler**



Rotor Tipleri;

Model		FC5000-1B	FC5000-4	T100-A/1	T100-A/2	T100-C
Rotor Tipi	Kapasite	RPM RCF(xg)	RPM RCF(xg)	RPM RCF(xg)	RPM RCF(xg)	RPM RCF(xg)
Sabit Açılı	12 x 0.5ml				16000 18320	16000 18320
Rotor	48 x 0.5ml					13500 18320
	12 x 1.5ml				16000 18320	16000 18320
	24 x 1.5ml					13500 18320
	4 x 5ml			10000 11000		
	12 x 5ml					15000 18320
	10 x 7ml					13500 18320
	4 x 10/15ml			10000 11000		
	24 x 10/15ml			10000 11000		
	6 x 50ml			10000 11000		
	4 x 100ml			10000 11000		
Açılır Rotor	24 x 5ml	4500 4360				
	40 x 5ml	4500 4360				
	24 x 7ml	4500 4360				
	60 x 7ml	4500 4360				
	120 x 7ml		3500 4850			
	88 x 5/7ml		3500 4850			
	8 x 10/15ml	5000 4360				
	12 x 10/15ml	5000 4360				
	16 x 10/15ml	5000 4360				
	24 x 10/15ml	5000 4360				
	28 x 10/15ml		5000 5180			
	32 x 10/15ml	5000 4360				
	64 x 10/15ml		4000 3400			
	72 x 10/15ml		3500 4850			
	32 x 15ml		3500 3400			
	36 x 20ml		3500 4850			
	8 x 50ml	5000 4360				
	16 x 50ml		3500 3400			
	4 x 100ml	5000 4360				
	8 x 100ml		4000 3400			
	4 x 500ml		3500 4850			
Microplate Rotor	2 x 2 x 96 x 0.2ml	5000 4360				
	4 x 2 x 96 x 0.2ml		5000 5180			