

FBC Serisi Yer Tipi Dokunmatik Panel Soğutmalı Santrifüj Rotor Tipleri



Yeni Nesil Dokunmatik Ekran



Blok Tipleri;

Model Rotor Tipleri	Kapasite (adet x ml)	FBC6000T RPM RCF(xg)	FBC-7000T RPM RCF(xg)	FBC-21000T RPM
RCF(xg)				
Sabit Açılı Rotor	12x10			21000 44440
	16x10			20000 46140
	8x20/30			20000 42040
	6x50			20000 44700
	8x50			15000 27540
	6x70			15000 25480
	8x100			12000 19610
	4x300	6000 5390		11000 18120
	6x300	6000 5660		10000 15730
	6x500	6000 6600		8000 11740
Açılır Rotor	4x5			15000 21670
	4x30			13000 23580
	4x50			10000 16240
	4x500			4000 3380
	6x1000	4200 5100	7000 11900	
	6x1100	4200 5100		
	6x2400		4000 5600	
	2x6x1000		4200 5600	
Plate Rotor	4x4x96	4000 2840		4000 3380
Sürekli Santrifüj Rotor	1000ml			17000 31390
Rotor	1800ml			14000 21710
	1800ml			8000 8600
	3000ml			8000 9450
Aralıklı Rotor	1000ml			17000 31390
	1200ml			12000 19350
	3000ml			10000 14890
Dikey Rotor	16x5			20000 38340
	8x30			18000 32650
Kan Torbası Adaptör	2x8x500ml		4000 5600	4000 5600

