


FBC Series Floor Type Touch Panel Refrigerated Centrifuge Rotor Types



New Technology Touch Screen



Block Types;

Model	Capacity	FBC6000T		FBC-7000T		FBC-21000T	
Rotor Type	(qty x ml)	RPM	RCF(xg)	RPM	RCF(xg)	RPM	RCF(xg)
Angle Rotor	12x10					21000	44440
	16x10					20000	46140
	8x20/30					20000	42040
	6x50					20000	44700
	8x50					15000	27540
	6x70					15000	25480
	8x100					12000	19610
	4x300	6000	5390			11000	18120
	6x300	6000	5660			10000	15730
	6x500	6000	6600			8000	11740
Swing Rotor	4x5					15000	21670
	4x30					13000	23580
	4x50					10000	16240
	4x500					4000	3380
	6x1000	4200	5100	7000	11900		
	6x1100	4200	5100				
	6x2400			4000	5600		
	2x6x1000			4200	5600		
Plate Bucket	4x4x96	4000	2840			4000	3380
Continuous Centrifugal Rotor	1000ml					17000	31390
	1800ml					14000	21710
	1800ml					8000	8600
Interval Rotor	3000ml					8000	9450
	1000ml					17000	31390
	1200ml					12000	19350
Vertical Rotor	3000ml					10000	14890
	16x5					20000	38340
	8x30					18000	32650
Bloodbag Adaptor	2x8x500ml			4000	5600	4000	5600

