


INORapid Test CORONAVIRUS (COVID-19)



INTERPRETATION OF RESULTS



As a supplementary detection

- To prevent missed detection
- Detect IgM and IgG antibodies of the coronavirus, and cooperate with nucleic acid detection to improve the detection rate of patients
- Easy sample collection ordinary labs can carry out inspections

Colloidal Gold Method

- Normally 15 mins, fastest 150s to get result
- Detect IgM and IgG, improve detection rate
- Supporting serum, and plasma sample types
- Optimizing N and S protein fragments, effectively improving the detection rate

Assay Kit by Colloidal Gold Method

Process	Specification	Product Code	Sample Loading
Colloidal Gold (IgM/ IgG) Method COVID-19 Antikor test kit	20 Test/Kit	INOPC 2200	Serum / Plasma: 10uL Whole Blood: 20uL